STARTERSGIDS

ROOD LICHTTHERAPIE

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- [1] Avci P, Gupta GK, Clark J, Wikonkal N, Hamblin MR. Low-level laser (light) therapy (LLLT) for treatment of hair loss. *Lasers Surg Med.* 2014 Feb;46(2):144–51. doi: 10.1002/lsm.22170. Epub 2013 Aug 23. PMID: 23970445; PMCID: PMC3944668.
- [2] Wunsch A and Matuschka K. A Controlled Trial to Determine the Efficacy of Red and Near-Infrared Light Treatment in Patient Satisfaction, Reduction of Fine Lines, Wrinkles, Skin Roughness, and Intradermal Collagen Density Increase. *Photomedicine and Laser Surgery.* Feb 2014; 32(2): 93–100.
- [3] Sadick NS. A study to determine the efficacy of a novel handheld light-emitting diode device in the treatment of photoaged skin. *Journal of Cosmetic Dermatology.* 2008 Dec.
- [4] Kim HK, Choi JH. Effects of radiofrequency, electroacupuncture, and low-level laser therapy on the wrinkles and moisture content of the forehead, eyes, and cheek. *Journal of Physical Therapy Science.* 2017 February.
- [5] *J Cosmet Dermatol.* 2007 Sep;6(3):189–94
- [6] Fan L, Yin R, Lan T, Hamblin MR. Photodynamic therapy for rosacea in Chinese patients. *Photodiagnosis Photodyn Ther.* 2018 Dec;24:82–87. doi: 10.1016/j.pdpdt.2018.08.005. Epub 2018 Aug 15. PMID: 30118905; PMCID: PMC6289730.
- [7] Zane C, Capezzer R, et al. Non-invasive diagnostic evaluation of phototherapeutic effects of red light phototherapy of acne vulgaris. *Photodermatology, Photoimmunology, and Photomedicine.* 2008 Oct.
- [8] da Silva Melo M, Alves LP, Fernandes AB, Carvalho HC, de Lima CJ, Munin E, Gomes MF, Salgado MAC, Zângaro RA. LED phototherapy in full-thickness burns induced by CO₂ laser in rats skin. *Lasers Med Sci.* 2018 Sep;33(7):1537–1547. doi: 10.1007/s10103-018-2515-8. Epub 2018 Apr 27. PMID: 29704068.
- [9] Baez F, Reilly LR. The use of light-emitting diode therapy in the treatment of photoaged skin. *Journal of Cosmetic Dermatology.* 2007 Sep.
- [10] Silveira PC, Ferreira KB, et al. Effect of Low-Power Laser (LPL) and Light-Emitting Diode (LED) on Inflammatory Response in Burn Wound Healing. *Inflammation.* 2016 Aug.
- [11] Alsharnoubi J, Shoukry K, et al. Evaluation of scars in children after treatment with low-level laser. *Lasers in Medical Science.*
- [12] Combination 830-nm and 633-nm Light-Emitting Diode Phototherapy Shows Promise in the Treatment of Recalcitrant Psoriasis: Preliminary Findings
- [13] Morita H, Kohno J. et al. Clinical Application of GaAlAs 830 nm Diode Laser for Atopic Dermatitis. *Laser Therapy.* 1993 vol 5.
- [14] Servetto N., Cremonuzzi D., Simes JC., Moya M., Soriano F., Palma JA., Campana VR. “Evaluation of inflammatory biomarkers associated with oxidative stress and histological assessment of low-level laser therapy in experimental myopathy”. *Lasers Surg Med.* 2010 August ;42(6):577–83. doi: 10.1002/lsm.20910.
- [15] Langella L., Casalechi H., Tomazoni S., Johnson D., Albertini R., Pallotta R., Marcos R., de Carvalho P., Leal-Junior E., “Photobiomodulation therapy (PBMT) on acute pain and inflammation in patients who underwent total hip arthroplasty—a randomized, triple-blind, placebo-controlled clinical trial”. *Lasers Med Sci.* 2018. Jun 16. doi: 10.1007/s10103-018-2558-x.
- [16] Lee JH., Chiang MH., Chen PH., Ho ML., Lee HE., Wang Yh. “Anti-inflammatory effects of low-level laser therapy on human periodontal ligament cells: in vitro study”. *Lasers Med Sci.* 2018; 33(3): 469–477. 2017 Nov 7. Doi: 10.1007/s10103-017-2376-6
- [17] Angelova A, Ilieva EM, et al. Effectiveness of High Intensity Laser Therapy for Reduction of Pain in Knee Osteoarthritis. *Pain Research and Management.* 2016;2016:9163618.
- [18] Fukuda VO, et al. Short-Term Efficacy of Low-Level Laser Therapy in Patients with Knee Osteoarthritis: A Randomized Placebo-Controlled, Double-Blind Clinical Trial. 2015 Dec 6;46(5):526–33.
- [19] Alayat MS, Aly TH, et al. Efficacy of pulsed Nd:YAG laser in the treatment of patients with knee osteoarthritis: a randomized controlled trial. *Lasers in Medical Science.* 2017 Apr;32(3):503–511.
- [20] Alfredo PP, Bjordal JM, et al. Efficacy of low level laser therapy associated with exercises in knee osteoarthritis: a randomized double-blind study. *Clinical Rehabilitation.* Jun 2012; 26(6): 523–33. at 635 nm for non-invasive body contouring of the waist, hips, and thighs. *Lasers in Surgery and Medicine.* Jan 2013;45(1):1–7.

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- [21] Effect of low-level laser therapy on muscle adaptation to knee extensor eccentric training. Baroni BM et al. 2015
- [22] Effects of Light-Emitting Diode Therapy on Muscle Hypertrophy, Gene Expression, Performance, Damage, and Delayed-Onset Muscle Soreness: Case-control Study with a Pair of Identical Twins. Ferraresi C et al. 2016
- [23] Low-level laser (light) therapy (LLLT) on muscle tissue: performance, fatigue and repair benefited by the power of light. Cleber Ferraresi et al. 2012
- [24] Does Phototherapy Enhance Skeletal Muscle Contractile Function and Postexercise Recovery? A Systematic Review. Paul A Borsa et al. 2013
- [25] Rossato M, et al. Time Response of Photobiomodulation Therapy on Muscular Fatigue in Humans. *Journal of Strength and Conditioning Research*. 2018 Nov;32(11):3285-3293.
- [26] Hofling DB, Chavantes MC, et al. Low-level laser in the treatment of patients with hypothyroidism induced by chronic autoimmune thyroiditis: a randomized, placebo-controlled clinical trial. *Lasers in Surgery and Medicine*. May 2013; 28(3): 743-53.
- [27] Zhao J, Tian Y, Nie J, Xu J, Liu D. Red light and the sleep quality and endurance performance of Chinese female basketball players. *J Athl Train*. 2012 Nov-Dec;47(6):673-8. doi: 10.4085/1062-6050-47.6.08. PMID: 23182016; PMCID: PMC3499892.
- [28] Figueiro MG, Plitnick B, and Rea MS. Light Modulates Leptin and Ghrelin in Sleep-Restricted Adults. *International Journal of Endocrinology*. 2012.
- [29] Paolillo FR, Borghi-Silva A, et al. New treatment of cellulite with infrared-LED illumination applied during high-intensity treadmill training. *J Cosmet Laser Ther*. 2011 Aug;13(4):166-71.
- [30] McRae E and Boris J. Independent evaluation of low-level laser therapy at 635 nm for non-invasive body contouring of the waist, hips, and thighs. *Lasers in Surgery and Medicine*. Jan 2013;45(1):1-7.
- [31] Sene-Fiorese M, Duarte FO, de Aquino Junior AE, Campos RM, Masquio DC, Tock L, de Oliveira Duarte AC, Dâmaso AR, Parizotto NA, Bagnato VS. The potential of phototherapy to reduce body fat, insulin resistance and "metabolic inflexibility" related to obesity in women undergoing weight loss treatment. *Lasers Surg Med*. 2015 Oct;47(8):634-42. doi: 10.1002/lsm.22395. Epub 2015 Jul 29. PMID: 26220050.
- [32] Núñez-Álvarez C, Del Olmo-Aguado S, Merayo-Llves J, Osborne NN. Near infra-red light attenuates corneal endothelial cell dysfunction in situ and in vitro. *Exp Eye Res*. 2017;161:106-115. doi:10.1016/j.exer.2017.06.008
- [33] Jin-Chul Ahn, Young-Hoon Kim, Chung-Ku Rhee. The effects of low level laser therapy (LLLT) on the testis in elevating serum testosterone levels in rats. *Biomedical Research* 2013; 24 (1): 28-32
- [34] European College of Neuropsychopharmacology (ECNP). "Lack of interest in sex successfully treated by exposure to bright light." *ScienceDaily*. ScienceDaily, 18 September 2016.
- [35] Liu KH, Liu D, et al. "Comparative effectiveness of low-level laser therapy for adult androgenic alopecia: a system review and meta-analysis of randomized controlled trials." *Lasers in Medical Science*. 2019 Aug.

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